



**oolab**  
ORTHOTICS

*Driving the orthotic revolution.*

# OOLAB PRESCRIPTION FORM

42 Niagara Street  
Hamilton, ON L8L 6A2  
Tel: 905-521-1230  
1-888-873-3316  
Fax: 905-521-1210

## DR/CLINIC INFORMATION

Name: \_\_\_\_\_  
Address: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
Account#: \_\_\_\_\_

## PATIENT INFORMATION

PLEASE PRINT  
First Name: \_\_\_\_\_  
Last Name: \_\_\_\_\_  
 M  F  Youth  
Weight: \_\_\_\_\_ (lbs)  
Shoe Size: \_\_\_\_\_ Age: \_\_\_\_\_

## PRESCRIPTION ORDER FORM

S/O # \_\_\_\_\_  
Date \_\_\_\_\_  
 Rush (\$25.00 fee)

### STANDARD ORTHOTIC

1mm  2mm  3mm  Poly

### DRESS ORTHOTIC

Carbon flex shell with shallow heel cup  
 1mm  2mm  Poly  
 Narrow Grind  Fashion Grind

### SOFT SPORT ORTHOTIC

Carbon flex shell, poron arch fill, 1mm underlay, & 1/8" ETC top cover  
 1mm  2mm  3mm

### PODIATRIC ORTHOTIC

3mm Polypropylene shell, soft arch fill, 1mm underlay, & 1/8" microcell top cover

### UCBL

Med/lat flanges, extrinsic heel posts, deep heel cup & microcell top cover  
 2mm  3mm

### CONTROL SPORT ORTHOTIC

(Add \$10.00)

Shell with dual density rearfoot post, skid plate, & 1/8" ETC top cover  
 2mm  3mm  Poly

### SOFT ORTHOTIC

(Add \$10.00)

Trilaminar shell material, 1/16" nylax forefoot extension, & 1/16" ultra hyde top cover

### GAIT PLATE

(Add \$10.00)

Poly shell, med/lat flanges & deep heel cup

Induce out toeing  Induce in toeing

### DIABETIC/ARTHRITIC

(Add \$10.00)

Carbon shell, 1/4" poron arch fill, 1mm underlay, & diabetic top cover  
 1mm  2mm

### XT CARBON

(Add \$35.00)

1.6mm Flex  2mm Semi-Flex  2.3mm Rigid

### TOP COVER

(Default Full)

3/4  Sulcus  Full

### TOP COVER MATERIALS

(Default ETC)

Black nylax can be substituted for blue poron

Black vinyl -->  Only  1/16"  1/8"

Blue vinyl -->  Only  1/16"  1/8"

Ultra Hyde -->  Only  1/16"  1/8"

1/8" ETC  Blue  Black

1/8" Microcell  Black  Grn/Blue  
 Pink/Purple

1/16" Suede  Tan  Brown  Black

Suede Only  Tan  Brown  Black

Plastazote 1/8"  1/16"

1/8" Neoprene  1mm Black Puff

### FOREFOOT CUSHION

1/16" nylax  1/8" nylax

### EXTRA CHARGE 1/8" TOP COVERS

(Add \$10.00)

Perforated Leather

Black Spenco

X-Static

Diabetic

### UNDERLAY MATERIALS

(Default Puff)

Full Foot  Forefoot  None

Vinyl  1mm Black Puff

Grey Suede  Durasole

### SHELL MODIFICATIONS

Heel Cup Height

Flat  Shallow  Standard  Deep

1st Ray Cutout  L  R

1st Met Cutout  L  R

Medial Flange (Shell)  L  R

Medial Flange (Soft)  L  R

Lateral Flange (Shell)  L  R

Lateral Flange (Soft)  L  R

Heel Hole  L  R

Rigid Morton's Ext (Poly only, Add \$10.00)  
 L  R

### ORTHOTIC ADDITIONS

Metatarsal Pad  L  R  
Move Distally \_\_\_\_\_ mm  Low Profile

Metatarsal Accom Pad  L  R  
Cutout L: 1  2  3  4  5   
Cutout R: 1  2  3  4  5

Metbar pad  L  R

Neuroma Padding Location:  2/3  3/4  
Sweet Spot Location: \_\_\_\_\_

Morton's Ext  Left  Right  
 Soft  EVA

Morton's Rev. Ext.  Left  Right

Heel Cushioning  Left  Right

Heel Spur Pad  Left  Right

Arch Fill  L  R  
 Soft  Firm  Hard  
 Low Profile  To Floor

Kinetic Wedge  Left  Right

Arch Cookie  1mm  2mm  3 mm

Rubber Skid Plate  Left  Right

### EXTRINSIC CORRECTIONS

Neutral  L  R

Heel Lift  L \_\_\_\_\_ mm  R \_\_\_\_\_ mm

Medial Wedge  L \_\_\_\_\_  R \_\_\_\_\_

Lateral Wedge  L \_\_\_\_\_  R \_\_\_\_\_

Varus  L \_\_\_\_\_  R \_\_\_\_\_

Valgus  L \_\_\_\_\_  R \_\_\_\_\_

Forefoot Medial  L \_\_\_\_\_  R \_\_\_\_\_

Forefoot Lateral  L \_\_\_\_\_  R \_\_\_\_\_

Extended Medial  L  R

### INTRINSIC CORRECTIONS

Varus  L \_\_\_\_\_  R \_\_\_\_\_

Valgus  L \_\_\_\_\_  R \_\_\_\_\_

### CAST DRESSING:

Min  Max

## SHOE ORDER

**1** SHOE: \_\_\_\_\_  
SIZE/WIDTH: \_\_\_\_\_  
COLOUR: \_\_\_\_\_  With Orthotic

**2** SHOE: \_\_\_\_\_  
SIZE/WIDTH: \_\_\_\_\_  
COLOUR: \_\_\_\_\_  With Orthotic

**3** SHOE: \_\_\_\_\_  
SIZE/WIDTH: \_\_\_\_\_  
COLOUR: \_\_\_\_\_  With Orthotic

### CLINICAL SUPPLIES

Brochures-Free

### IMPRESSIONS

RX Forms-Free

6-Pack-\$38.00

\*Shipping and HST

25-Pack-\$125.00

Extra Charge

50-Pack-\$199.00

50-Plaster Splints-\$58.00

### ADDITIONAL INSTRUCTIONS

#PAIR ORTHOTICS ONLY \_\_\_\_\_

#PAIR SHOES ONLY \_\_\_\_\_

#PACKAGE ONLY \_\_\_\_\_

(1 Pair Orthotics/1 Pair Shoes)

